

OPERATING INSTRUCTIONS

Hookup

Connect the Probe to the unit with the fifteen-foot (15') Probe Cord. Lubricate the Probe; elevate the tail until it is horizontal, not up over the back. Press the Probe into the rectum with a slow positive push. Aim the nose of the Probe downward to facilitate entrance, aligning the path of the Probe with the slant of the floor of the pelvis. The Probe electrodes on the two-electrode Probe should be ventrally oriented and the handles pointing upward.

Operation

Be sure that the Power Control (left) is turned all the way counter-clockwise, then put the Function Switch in the BATT position. From this point through a normal ejaculation, only the Power Control is used. Efficient and consistent results are dependent upon development of a satisfactory technique of manipulating the Power Control. This is best accomplished by careful observation of the animal's reaction and proper correlation of the relative meter indication in order to place the frequency of the stimulus and to ascertain the power level required for each particular animal being stimulated. Apply the stimuli, rhythmically, according to how the animal is reacting.



CAUTION: If you are using one of our new Probes as pictured in this manual, be careful about turning the Power Control up too far at first. Proceed very cautiously as you can expect to get a very prompt protrusion and erection followed by an immediate ejaculation, all possibly in 30 seconds (according to many reports from your colleagues who got one of the first Thousand of the new Probes).



Observe the underside of the belly; watch for movement indicating position of the penis. Apply stimuli, up to a given power setting, for 5 to 10 impulses (one-half minute or so) if there is no evidence of enlargement of the penis and forward movement through the sheath during this cycling, increase the power level slightly and proceed through another series of 5 to 10 impulses. If penis has not yet become visible, or it is just protruding from among the preputial hairs, reach up behind and above the neck of the scrotum with your strongest arm. (Use the bull's leg, above the hock, as a fulcrum for your arm.) Turn the fingers into the hand so they are bent at right angles of the first joint, and force these fingers into the neck of the scrotum, just like the action of a rake! Work your finger rake into the tissue until you feel the curvature of the sigmoid flexure, then force your hand forward until you feel the flexure uncoil. All the time you are performing this feat, it will help if you can maintain the rhythmic stimulations. Once the sigmoid flexure is straightened out, and the penis visibly extruded and beginning to take on a pink color, it is a matter of keeping the penis extruded and erect for a few more moments while increasing the power of the stimuli steadily until ejaculation occurs.

After ejaculation occurs, rotate the Power Control to its counter-clockwise position and place the Function Switch in the OFF position .